



## Compressor Nebulizer

Ideal choice for achieving patients compliance



The MiniCare is the brand new portable compressor nebulizer just introduced to the market with the best interest of yours in mind so you can be at ease knowing you are in good hands with proper care of your aerosol medication therapy. It is engineered for durability and cost effectiveness. The efficient nebulizer is just priced right with high quality performance when cost is a consideration in your mind. It is an ideal alternative to replace other expensive nebulizer designs.

**APEX**

## Mini-Care



**Respiratory Therapy**

Care for a Healthy Life

### Key Features:

- ◆ Economically priced to meet the needs of the homecare market.
- ◆ Efficient aerosol delivery; shorter treatment time.
- ◆ Attractive modern dedicated compartment space keep the medications and accessories safe, ensuring all the necessities are at hand when required.
- ◆ Piston-driven compressor assures strong steady pressures and flows.
- ◆ Unique nebulizer kit holding slot is built into the unit for ease operation.
- ◆ Quiet, low-noise performance provides unobtrusive aerosol therapy.
- ◆ The power cord tucks neatly into the back of the unit.
- ◆ The sturdy built-in handgrip makes it easy to carry.
- ◆ Durable, reliable and solid state components.



Nebulizer Holder



Large Storage Compartment

Hidden Storage for Electrical  
Cord

Convenient built-in Handgrip

### Specifications

Specification	Description
Electrical Requirement:	AC 100/120V, 60 Hz, 1.2A ( for 120V system) AC 230V 50/60 Hz, 0.6A ( for 230V system)
Maximum Pressure:	≅ 30 psi
Maximum Flow Rate:	≅ 12 lpm
Operating Flow Rate:	≅ 5 lpm
Nebulization Rate:	0.2-0.3ml/min
Particle Size	0.5-5 microns*
Dimensions (W x D x H):	29.8 x 17.3 x 10.7 cm / 11.7" x 6.8" x 4.2"
Weight:	2.1 kg / 4.6 lb

\*Particle size and nebulization rate are determined by the nebulizer used with the device


All specifications are subject to change without notice.



## APEX MEDICAL CORP.

9, Min Sheng St, Tu-Cheng,  
Taipei County, 236, Taiwan, R.O.C.  
TEL: 886-2-22683100  
FAX: 886-2-22686525  
E-mail: marketing@apexmedical.com.tw  
http://www.apexbrand.com

Distributed by:

 <b>IHS</b> Independent Healthcare Service Ltd. Unit 3 Westpoint Business Park, Mulhuddart, Dublin 15.	<b>SALES &amp; SERVICE</b>  <b>Tel: 01 443 3390</b>  <b>Email: info@ihsltd.ie</b>  <b>Website: www.ihsltd.ie</b>  <i>Working Together With Care...</i>
---	--