

Eva600EE™

SystemRoMedic



Lifting capacity 270 kg / 600 lbs • Extra strong and safe

The Eva600EE is a versatile mobile patient lift for all lifting situations, especially for heavier patients. Together with SystemRoMedic's functional and easy-to-use accessories, the Eva600EE is the right choice for lifting in connection with transfers to and from the floor, chair or bed, as well as for horizontal lifting and gait training.

Strong and safe for heavy patients:

- Lifting capacity 270 kg / 600 lbs
- Aluminium, for light weight
- Adjustable handle
- Emergency stop
- Electrical and manual emergency lowering
- Battery with long service life
- Built-in charging
- Mast adjustable to three height positions
- Maintenance-free wheels
- Electrical base-width adjustment
- Wide slingbar
- Many accessories



Eva600EE
complies with requirements for Class 1-products in the European Medical Device Directive MDD/92/43/EEC.

Product specifications:

Eva600EE _____ Art. No. **60100003**

TECHNICAL INFORMATION

| | |
|------------------------------------|--------------------------|
| Safe working load: | 270 kg / 600 lbs |
| Lifting height min-max: | 500-1810 mm (19,7-71,2") |
| Outer dimensions base: | 750-1150 mm (29,5-45,2") |
| Inner dimensions base: | 650-1050 mm (25,6-41,3") |
| Lifting speed <i>without</i> load: | 3.5 cm/s (1,4"/s) |
| Emergency lowering: | Manual and electrical |
| Castors: | Ø 100 mm (Ø 3,9") |
| Battery capacity: | 2.9 Ah |
| Material: | Aluminium |
| Protection class: | IPx4 |
| Sling bar (width): | 600 mm (23,6") |
| Weight (with battery): | 39 kg / 85,9 lbs |
| Weight (without battery): | 37 kg / 81,5 lbs |



The Eva600EE has a sturdy design and is made entirely of aluminium. The lift is strong and stable, light-weight, easy to manoeuvre and easy to keep clean. Aluminium does not rust and the material is also environmentally friendly, as it can be recycled multiple times.



The Eva600EE has an emergency stop button easily accessible on the control box. The lift also has both electrical and manual emergency lowering.

The battery has high capacity and will perform many lifts per charge, which means less frequent charging. An indicator on the control box shows when charging is in progress and when the battery is fully charged. Built-in charging is standard on Eva600EE.

The mast's three different height positions give the Eva600EE a very large lifting range. The lift can be easily adapted for different lifting situations, for example, for lifting from the floor and for gait training.

The ergonomic handle is height-adjustable to enable the caregiver to work more conveniently, comfortably and in an ergonomically correct manner that prevents occupational injuries and minimizes discomfort for those who already have an injury. The handle's smooth surface is easy to keep clean.



Encapsulated dual wheels.

The Eva600EE has encapsulated dual wheels which provide for a smoother ride for the user while making the lift very easy to roll, even on softer surfaces. The encapsulated wheels do not gather dust and dirt from the floor, making them low-maintenance.

The hand control has raised, easy-to-push buttons for manoeuvring the lift and adjusting the base width.

The Eva600EE comes with a wider sling bar (600 mm), providing a more open sling which is ideal for larger users.



SystemRoMedic

"We love easy transfers"



The philosophy behind SystemRoMedic is simple. We love to develop and manufacture assistive devices that make life easier for both users and personnel in the care sector.

SystemRoMedic comprises four product areas that meets most needs in different transfer situations.

- Transfer**, between two locations.
- Positioning**, within the same location.
- Support**, for mobility.
- Lifting**, both manual and mechanical.

SALES & SERVICE

Independent Healthcare Service Ltd.
Unit 3 Westpoint Business Park,
Mulluddara, Dublin 15.

Tel: 01 443 3390

Email: info@ihsltd.ie

Website: www.ihsltd.ie

Working Together With Care...