



Enjoy a great day on its watch!



Fingertip SpO<sub>2</sub> incorporates the electronics and sensor into one unit that provides cost-effective solution for spot-checks and short-term monitoring.

Fingertip SpO<sub>2</sub> is not influenced by patient's motion; therefore, it gives patients the freedom of having their physical checkup taken anywhere and anytime they want including having it incorporated into an athletic activity.



**Respiratory Therapy**

**APEX**

**SB 100**



Carefor a Healthy Life

### Features/Benefits:

- ◆ Delivers accurate pulse rate and blood saturation in seconds
- ◆ One-touch keypad for easy operation
- ◆ Big and bright LED display
- ◆ Compact and light for mobility
- ◆ Two AAA Alkaline batteries for easy power supply
- ◆ Automatic power-off after 10 seconds in idling



### Technical Specification

Scope of Measurement	SpO2 : 35% ~ 99% Pulse rate: 30 ~250 bpm
Precision	SpO2 : 70% ~ 100% ± 2 digits(%) 35% ~ 69% : unspecified
Pulse rate	30 ~ 250±3 bpm
Averaging	4 for SpO2 ; 8 for Pulse rate
Weight	37g
Size(HxWxD)	63.5 x 34 x 36 mm
Power Supply	2 "AAA" alkaline (Can be used for 14 hours)
Environmental Conditions	Operating Temperature: 41°F ~ 104°F (5°C ~ 40°C) Storage Temperature : -4°F ~ 158°F (-20°C ~ 70°C) Operating Humidity : 15% ~ 95%, noncondensing
Standards	Type of Protection : IEC60601-1-2, Class B , IEC60601-1, Type BF , ISO 9919: 2005

\* All specifications are subject to change without notice.

## APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,  
Taipei County, 236, Taiwan, R.O.C.  
TEL: 886-2-22683100  
FAX: 886-2-22686525  
E-mail:marketing@apexmedical.com.tw  
http://www.apexbrand.com

Distributed by:

	<b>SALES &amp; SERVICE</b>
<small>Independent Healthcare Service Ltd.</small>	<b>Tel: 01 443 3390</b>
<small>Unit 3 Westpoint Business Park, Mulhuddart, Dublin 15.</small>	<b>Email: info@ihsltd.ie</b>
	<b>Website: www.ihsltd.ie</b>
<i>Working Together With Care...</i>	