



## Health at fingertip



If a handheld device at a very attractive price for monitoring patient's oxygen saturation in blood and heart rate with precision in all ranges from pediatric to adult, please look no further, Apex is the one.



**Respiratory Therapy**

Care for a Healthy Life

### Features/Benefits:

- ◆ Touch-keypad for easy operation
- ◆ Backlight display for night or day
- ◆ High & low SpO<sub>2</sub> and pulse rate programming with audible/visual alarm
- ◆ Sound volume adjustment
- ◆ Pock-size design
- ◆ 4 AA Alkaline batteries for power supply
- ◆ Auto-off to conserve battery life
- ◆ Data management software ; 3-mode versatility: Spot checking, memory, and PC-Link
- ◆ Orange pouch for the protection and mobility of the oximeter and probe



### Technical Specification

Scope of Measurement	SpO <sub>2</sub> : 35% to 100% Pulse rate: 30 - 250 bpm
Precision	SpO <sub>2</sub> : 70% ~ 100% ± 2 digits(%) 35% ~ 69% : unspecified
Pulse rate	30 ~ 250 ± 3 bpm
Data Update period	less than 2 seconds
Averaging	4 for SpO <sub>2</sub> ; 8 for Pulse rate
Weight	280g ; 0.33 ponds (excluding battery and probe)
Size(HxWxD)	14.5 x 7.25 x 2.25 cm /5.70" x 2.85" x 0.89"
Battery	4 "AA" alkaline /15 hrs continuous operation.
Environmental Conditions	Operating Temperature: 32°F ~ 113°F (0°C ~ 45°C) Storage Temperature : -4°F ~ 122°F (-20°C ~ 50°C) Operating Humidity : 15% ~ 90%, noncondensing
Standards	Type of Protection : Class II equipment / IEC60601-1-2, Class B / EN60601-1, Type BF / ISO 9919: 2005

\* All specifications are subject to change without notice.

## APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,  
Taipei County, 236, Taiwan, R.O.C.  
TEL: 886-2-22683100  
FAX: 886-2-22686525  
E-mail: marketing@apexmedical.com.tw  
http://www.apexbrand.com

Distributed by:

<b>IHS</b>	<b>SALES &amp; SERVICE</b>
<small>Independent Healthcare Service Ltd.</small>	<b>Tel: 01 443 3390</b>
<small>Unit 3 Westpoint Business Park, Mulhuddart, Dublin 15.</small>	<b>Email: info@ihsltd.ie</b>
	<b>Website: www.ihsltd.ie</b>
<i>Working Together With Care...</i>	