



Independent Healthcare Service Ltd.



Support Surfaces

PM170

5" Alternating Pressure Relief Mattress Overlay System

The PM170 is a high quality, affordable Alternating Pressure Relief Mattress Overlay system which is designed for pressure ulcer prevention and treatment used in basic care or moderate care situations. It provides alternating pressure therapy with low air loss features. The Static mode provides a consistent pressure in all air cells to facilitate nursing care, dressing changes and patient repositioning. The low air loss laser holes in the mattress provide adequate ventilation and moisture control. The mattress cover is waterproof, vapor permeable and low friction, reducing moisture and shear force. The pump operates quietly and without vibration to allow the patient to rest comfortably without interruption.



Support Surfaces



Features:

- High Performance air pump with 8 liter compressor.
- Quiet, Vibration-free pump operations.
- Convenient switch to change between alternating therapy and static mode.
- 8-minute cycle time (120V 60 Hz power supply).
- 9.6-minute cycle time (230V 50 Hz power supply).
- Conspicuous visual alerts.
- 17 structurally independent air cells facilitate repair, replacement and cleaning.
- Top 3 air cells remain inflated to serve as head pillow.
- CPR tags.
- Low air loss.
- Waterproof, vapor-permeable, low friction cover.
- Full periphery zipper on cover allows complete detachment from mattress.



APEX quick connector provides easy, secure connection between pump and mattress.



User friendly control panel with various comfort settings.



Easy to use hanger for attaching the pump to the bed frame.

Specifications:

PUMP	Power Supply	AC120V/60Hz 1.6A (for 120V system) AC230V/50Hz 0.06A
	Pressure Range (mmHg)	40-80 mmHg
	Cycle Time (min)	8 min / 60Hz, 9.6 min / 50Hz
	Dimensions (L x W x H)	29 x 18.5 x 12.6 (cm) or 11.4" x 7.3" x 5.0"
	Weight	2.2 Kg or 5 lb
MATTRESS	Arrangement	17 x 5" air cells
	Material	Cover Nylon/PU Cells Nylon/PVC Base Nylon/PVC
	Dimensions (L x W x H)	200 x 90 x 12.7 (cm) or 78.7" x 35.4" x 5"
	Weight	5.2 kg or 11.4 lb



Independent Healthcare Supplies

Distributed by:



Independent Healthcare Service Ltd.
Unit 3 Westpoint Business Park,
Mullhuddart, Dublin 15.

Tel: 01 443 3390
Email: info@ihsltd.ie
Website: www.ihsltd.ie

Working Together With Care...