

# Bedroom Equipment



A complete range  
for hospital rooms

A harmonious environment  
for your patients



# Hospital Room

The  
professional  
solution



The Hospital room has been designed for clinics and hospitals. Its equipment is both technical and functional, to meet the needs of both patients and hospital departments. A professional solution suited to purpose.

**Hospitec bed** (4952)



**Boléro Armchair and foot rest** (6120 001) (6189 021)



**Boléro chair** (6187 005)



**Boléro bridge** (6188 005)



**Lateral release table** (5647 000)



**Bedside cabinet** (5400)



# Bariatric range

Comfort and  
Safety for  
the patient



## Atlas bed

Lifting capacity: 320 kilos  
Height adjustable from 35 to 66 cm  
Width: 120 or 140 cm  
Length: 200 cm



## Boléro bridge

(6198000)  
Height: 85 cm  
Width: 78 cm  
Depth: 50 cm  
Seat height: 52 cm  
Weight capacity: 270 kg



**Boléro  
armchair** (6128000)

**Weight capacity: 270 kg**

Height: 118 cm

Width: 79 cm

Depth: 52 cm

Seat height: 53 cm

Height of removable/adjustable  
armrests: 70 cm

**Variation :**

Height of fixed armrests: 72 cm

**Bedside cabinet**

(5500000)

Height: 83 cm

Width: 48 cm

Depth: 51 cm



**Table**

(5648000)

Height: 78 cm

Width: 110 cm

Depth: 60 cm

# Bedside cabinets



## Shared features (5300 series)

- Cabinet on 4 castors
- Low-level bumper to protect walls
- Removable ABS cover RAL GREY 7035 protective top, with rim

### Medical bedside cabinet

version 1  
(5300 000)

Drawer with removable plastic container.  
2-door facade.

Height: 86 cm  
Width: 54 cm  
Depth: 56 cm



### Bedside cabinet with refrigerator

(5300 007)

Door and drawer facade  
Refrigerator capacity:  
30 litres.

Height: 86 cm  
Width: 54 cm  
Depth: 56 cm



### Medical bedside cabinet version 3 (5300 010/030)

Height: 86 cm  
Width: 54 cm  
Depth: 56 cm



Door (right- or left-opening).  
Drawer with removable plastic container on facade.  
Sliding telephone shelf.  
Glass holder and bottle holder on side facade.  
ABS container with basin holder and shelf integrated into the back.



Optional extra:  
removable  
chrome rail

## Medical bedside cabinet (5400)

Height: 87 cm  
Width: 48 cm  
Depth: 45 cm



### Version 1

Sliding cabinet with storage at 2 levels.  
Space phone under upper tray.  
Rail along 2 side of top surface.  
Basin holder stores vertically in the side.



### Version 2

1 drawer, 1 door.  
Door opening to left or right (to be specified at time of ordering).  
Storage area under the top surface.  
Rail along 1 side of top surface.  
Choice of 2 castors / 2 feet or 4 castors.



Height: 79 cm  
Width: 42 cm  
Depth: 46 cm

# Seating



**Ergonomic Design**

**Comfort for your patients  
and users**

**Technological sophistication**



# Boléro Range

Comfort and Hygiene

## Technical specifications

- Tubular Structure, epoxy finish: RAL GREY 7035
- PVC M1 coated, Griffine, 'Diabolo Club' Range
- One piece back section (standard version ), option 2/3 - 1/3 and optional removable height-adjustable headrest
- Bilateral fitments for IV drip rod and urine bag holder

### Armchair (6120 001)

Height: 108 cm  
 Width: 66 cm  
 Seat width: 50 cm  
 Seat depth: 49 cm  
 Depth: 47 cm  
 Seat height: 50 cm  
 Armrest height: 74 cm



### Foot rest (6189 021)

Adjustable incline  
 Height: 42 cm  
 Width: 51 cm  
 Depth: 35.5 cm  
 Seat height: 35 cm

**Variation** : fixed foot rest

### Chair (6187 005)

Height: 85 cm  
 Width: 46 cm  
 Depth: 60 cm  
 Seat height: 49.5 cm  
 Seat depth: 40 cm



### Bridge (6188 005)

Height: 85 cm  
 Width: 58 cm  
 Depth: 60 cm  
 Seat height: 49.5 cm  
 Seat depth: 40 cm



## Reclining Easychair (6120 004)



- Tubular Structure, epoxy finish: RAL GREY 7035
- M1 PVC upholstery
- 4 feet
- Incorporated headrest
- Optional 2/3 - 1/3 backrest
- Optional removable/height-adjustable headrest
- Optional accessories: depth-adjustable tray, IV drip rod, leg-rest, urine bag holder
- Seat height: 47 cm
- Seat depth: 49 cm
- Width between armrests: 52 cm
- Bilateral fitments for IV drip rod and urine bag holder

## Synchrone Easychair (6120 005)

- Tubular Structure, epoxy finish: RAL GREY 7035
- M1 PVC coated
- 4 patins
- Incorporated headrest
- Optional 2/3 - 1/3 backrest
- Optional removable/height-adjustable headrest
- Optional accessories: depth-adjustable tray, IV drip rod, leg-rest, urine bag holder
- Seat height: 47 cm
- Seat depth: 49 cm
- Width between armrests: 52 cm
- Bilateral fitments for IV drip rod and urine bag holder

