

Pediatrics




Bambino Bed

ref. 4700 L S1 2 PE

The 4 to
10 years old
children bed



Folding "Plexiglas" side rails
horizontal fitting (epoxy finish),



4 braked castors \varnothing 100 mm

Two-coloured melamineted
base structure

Model shown: Mistigri H/F board, base structure & side rails "Plexiglas"



Mattress platform:

- Metallic slats

Technical specifications

- Electric equipment with ingress protection : IP66
- Electric variable height (crosspiece system), 6 000 Newton:
 - Minimum height: 29.5 cm
 - Maximum height: 89.5 cm,
- Tubular base section, (50 x 30) and (40 x 30 mm), epoxy finish: RAL GREY 7035
- Adjustable wall saver,
- 4 braked castors Ø 100 mm
- Tubular platform, 40 x 30 mm, mattress 160 x 80 cm, epoxy finish. It is equipped with:
 - Welded metal slats
 - Bilateral fitments for lifting pole,
- Electric backrest, 70° angle , 2 500 Newton,
- Manual leg elevation by double ratchets, angle 20°

Options and Accessories

- Lifting pole (Epoxy finish),
- Fitting for IV Rod with adapter on fitting for lifting pole,
- Folding metallic side rails horizontal fitting (epoxy finish), complying with the NF EN 747-1 standards (Free space between the rungs from 60 to 75 mm),
- Folding "Plexiglas" side rails horizontal fitting (epoxy finish), complying with the NF EN 747-1 standards (full & transparent),
- Central locking for 4 castors, without directional castor,
- Electric leg elevation,
- Electric knee break,
- "Liberty" mattress, reference 4064 002,
- Two coloured "BAMBINO" or "MISTIGRI" head and foot boards,
- Two-coloured melaminated base structure

Crib

ref. 4 721 000

Newborn



Options and Accessories

- Crib Mattress

Cot

ref. 4 708 006

For age
 6 months
 to 4 years



Options and Accessories

- Cot Mattress
- Epoxy colour:



Soleil



Océan



Coquelicot