

Seating



Ergonomic Design

**Comfort for your patients
and users**

Technological sophistication



Boléro Range



Technical specifications

- Tubular Structure, epoxy finish: RAL GREY 7035
- PVC M1 coated, Griffine, 'Diabolo Club' Range
- One piece back section (standard version), option 2/3 - 1/3 and optional removable height-adjustable headrest
- Bilateral fitments for IV drip rod and urine bag holder

Armchair (6120 001)

Height: 108 cm
 Width: 66 cm
 Seat width: 50 cm
 Seat depth: 49 cm
 Depth: 47 cm
 Seat height: 50 cm
 Armrest height: 74 cm



Foot rest (6189 021)

Adjustable incline
 Height: 42 cm
 Width: 51 cm
 Depth: 35.5 cm
 Seat height: 35 cm
Variation : fixed foot rest

Chair (6187 005)

Height: 85 cm
 Width: 46 cm
 Depth: 60 cm
 Seat height: 49.5 cm
 Seat depth: 40 cm



Bridge (6188 005)

Height: 85 cm
 Width: 58 cm
 Depth: 60 cm
 Seat height: 49.5 cm
 Seat depth: 40 cm



Reclining Easychair (6120 004)



- Tubular Structure, epoxy finish: RAL GREY 7035
- M1 PVC upholstery
- 4 feet
- Incorporated headrest
- Optional 2/3 - 1/3 backrest
- Optional removable/height-adjustable headrest
- Optional accessories: depth-adjustable tray, IV drip rod, leg-rest, urine bag holder
- Seat height: 47 cm
- Seat depth: 49 cm
- Width between armrests: 52 cm
- Bilateral fitments for IV drip rod and urine bag holder

Synchrone Easychair (6120 005)

- Tubular Structure, epoxy finish: RAL GREY 7035
- M1 PVC coated
- 4 patins
- Incorporated headrest
- Optional 2/3 - 1/3 backrest
- Optional removable/height-adjustable headrest
- Optional accessories: depth-adjustable tray, IV drip rod, leg-rest, urine bag holder
- Seat height: 47 cm
- Seat depth: 49 cm
- Width between armrests: 52 cm
- Bilateral fitments for IV drip rod and urine bag holder



Ergo comfort

ref. 6 122 005

Mobility and
comfort



Model shown:
2/3-1/3 backrest
and centralised locking



Meal Tray



Removable / height
adjustable headrest



Anatomical reducer



Leg-rest

And also

- 4 braked castors 125 mm Ø
- Moulded seat
- M1 class anatomical upholstery, 'Griffine Diabolo Club' range
- Choice of colours (see colour chart)
- Space between the armrests: 53 cm
- Seat height: 52 cm
- Seat depth: 48 cm
- 2 push-handles
- Fitment for leg-rest
- Bilateral fitments for IV drip rod and urine bag holder

Options and Accessories

- Centralised locking
- Integral design castors if centralised locking
- 2/3-1/3 backrest
- Removable / height adjustable headrest
- Backrest 2/3-1/3 et tête amovible
- Synchronised seat and back section during reclining
- Anatomical reducer
- Leg-rest
- Integrated chest supports
- ABS Tray
- IV drip rod
- Urine bag holder

Navigation Chair

ref. 6 956 000

**Simplifying
Transfer**



Options and Accessories

- IV drip rod
- Document/File holder
- Smart lock system

Overall Dimensionst

Length: 93 cm
 Width: 68 cm
 Height: 112 cm

Seat

Width: 56 cm
 Depth: 43 cm
 Seat height: 55 cm

Backrest

Width: 49 cm
 Height: 40 cm



Side view



Smart lock system



Document/File holder



Centralised locking bar