



Independent Healthcare Supplies

new

Showering & Bathing

AquaMaster Tilt in Space Shower Chair

Chiltern Invadex has been helping to provide bathing solutions for domestic, healthcare and leisure environments for over 25 years.

The AquaMaster Tilt in Space Shower Chair is the newest addition to the range. It has been developed in conjunction with healthcare professionals to create a chair which provides comfort and security for the user combined with ease of use.

It provides a tilting angle of 0 to 30 degrees with an additional 10 degrees provided by the ergonomic angled backrest.

The chair boasts an electrostatically coated stainless steel frame which provides a durable, corrosion resistant finish and a dual gas spring assembly which allows smooth operation via a single rear mounted control handle.

It is available in two standard widths, 18" and 21" complete with an extensive range of seat, armrest, footrest and backrest padding configurations which provide a safe, user specific chair. This is further complemented by a complete range of accessories including fully adjustable head supports, foot trays, calf straps, side padding, harnesses and commode pans.

Seat, armrest and footrest configurations are retrospectively replaceable to allow the chair specification to be amended to suit degenerative conditions and to enable re-issue.

The AquaMaster Tilt in Space shower chair has been competitively priced, please contact Customer Services on 01869 365500 opt 1 for a quotation or to arrange a no obligation assessment.



Coming Soon
Tilting Shower Cradle



	Seat Width	Overall Width	Seat Height	Overall Height	Seat Depth	Overall Depth	Overall Tilted Depth	Safe Working Load
A14 18"	450	568	610	1380*	435	1065	1440*	160kg/25st
A15 21"	520	670	610	1380*	435	1065	1440*	

* complete with fabric head support at maximum extension





Independent Healthcare Supplies

Options and Accessories



Straight, half wide and wide armrests available as an option. Standard supplied with tip up armrests.



Seats available in a variety of aperture sizes, with extended options, skirts and bespoke seating available on request.



Full Support Harness



Padded Clip on Foot Tray



Commode Pan (round or oval) & Cradle



Gutter Armrest

Thoracic Support

Adjustable Head Support



This is a selection of options and accessories, many more are available.

Features and Benefits

- Designed and field tested in conjunction with Occupational Therapists and Healthcare Professionals.
- Offers a safe, controlled and comfortable showering position for those with severe physical conditions.
- Full range of aperture, horseshoe and flat padded seats available.
- Option for bespoke seating, including pressure relief foam.
- Armrest, footrest and backrest options.
- Extensive range of accessories available.
- Dual gas rams to ensure smooth reliable operation.
- Single hand control lever.
- Stainless steel, electrostatically powder coated construction for strength and durability.
- Corrosion resistant stainless steel and nylon braked castors. Directional locking option available.
- Commode cradle options for round or oval pan types.
- Also available as a tilting shower cradle.
- Full 12 month warranty.
- Free, no obligation technical survey.
- Carries the CE mark in accordance with the Medical Devices Directive.

For further details or to arrange an assessment visit please contact:

SALES & SERVICE

Tel: 01 443 3390

Email: info@ihsltd.ie

Website: www.ihsltd.ie

Working Together With Care...

AquaMaster Tilt in Space Shower Chair



BS EN ISO 9001:2008
FM553189