



Respiratory Therapy

Portable Suction Unit VacPlus

The VAC-series pump is Apex Medical's most recent generation of portable suction units and has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements. The VacPlus unit is lightweight, compact and can be used in either AC or DC mode. This unit is available with a rechargeable battery, remote controller.





Features:

Smooth hygienic finish, water resistant, easy to clean.



Integral secure canister holder.



Portable design with ergonomic easy-to-carry handle.



Pressure regulator control and angled, easy to read analogue pressure gauge.



Optional attached remote controller for additional operational convenience.

Accessories

Optional accessories available



Single-patient use canister available in 800ml or 1800ml



Filters and Tubing



Charger and remote controller

Specifications:

Electrical Rating	AC 100~240V,50/60Hz DC 12V
Power Consumption	36W
Pressure Range	0~600 mmHg (adjustable)
Maximum Flow Rate	24 LPM
Net Weight	3.3 kg
Dimension	352(L) x 206(W) x 192(H) mm.
Sound Level	< 55 dB
1 year warranty	



Approval:

- EN 60601-1
- EN 60601-1-2
- EN/ISO 10079-1



Independent Healthcare Service Ltd.

Distributed by:



SALES & SERVICE

Tel: 01 443 3390

Email: info@ihsltd.ie

Website: www.ihsltd.ie

Independent Healthcare Service Ltd.
Unit 3 Westpoint Business Park,
Mullhuddart, Dublin 15.

Working Together With Care...