



Respiratory Therapy

Portable Suction Unit VacPro

The VAC-series pump is Apex Medical's most recent generation of portable suction units and has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements. The VacPro is durable, lightweight and versatile in the provision of general suction usage.

Lightweight!
**1.8kg
only**





Features:

Smooth hygienic finish, water resistant and easy to clean.



Integral secure canister holder.



Lightweight and portable design with ergonomic easy-to-carry handle.



Pressure regulator control and angled, easy to read analogue pressure gauge.



Accessible and secure power connection.

Accessories

Optional accessories available



Single-patient use canister available in 800ml or 1800ml



Filters and Tubing.

Specifications:

Electrical Rating	AC 110~120V, 60Hz AC 220~240V, 50/60Hz
Power Consumption	90W
Pressure Range	0~600 mmHg (adjustable)
Maximum Flow Rate	24 LPM
Net Weight	1.8 kg
Dimension	352(L) x 206(W) x 192(H) mm.
Sound Level	< 55 dB
1 year warranty	



- Approval:
- EN 60601-1
 - EN 60601-1-2
 - EN/ISO 10079-1



Independent Healthcare Service Ltd.

Distributed by:



SALES & SERVICE

Independent Healthcare Service Ltd.

Unit 3 Westpoint Business Park,
Mulhuddart, Dublin 15.

Tel: 01 443 3390

Email: info@ihsltd.ie

Website: www.ihsltd.ie

Working Together With Care...