

The ideal aid for easy patient transport.



*Integrated patient scales (accessory)*

*The **GK hektor** chair-lifter with digital weighing machine can be used to weigh residents precisely and conveniently without their needing to change seats additionally.*

*Complete with TARE function and automatic switching off.*

*Toilet and bed-pan*

*The ergonomically shaped seat opening enables the **GK hektor** chair-lifter to be used as a toilet seat. A holder for a bed-pan under the seat is available as an optional extra.*



*Power supply designed for optimum use*

*The rechargeable battery can be detached for recharging at the external charger. With a reserve battery (accessory), the **GK hektor** chair-lifter is always ready for instant use.*



*Technical Specifications and Design subject to change without notice.*

# Hektor

## CHAIR-LIFTER



# Hektor

02/2005 - 1.000



OQS Reg. No. 1541/0  
EN ISO 9001



For the well-being of patient and nursing personel.



*The **GK hektor** chair-lifter – for bathing, showering, toilet, moving, laying down, weighing.*

*The narrow chassis and the broad, ergonomically designed convex hygienic seat of soft, skinfriendly polyurethane enable residents who can walk to sit down in comfort. Less mobile residents can be taken directly from their wheelchair or bed and can then be transported gently and safely to the bath-tub, the shower, or the WC.*

*The nursing staff are not exposed to health risks associated with strenuous lifting. The low chassis ensures that it can be driven under nursing beds. A detachable back support and all-round hand-holds provide the resident with additional comfort and security.*


*The **GK hektor** chair-lifter is ready for use immediately it is delivered.*



*The perfect extra for use with GK lifting bathtubs.*



*Height adjustable from floor level to 600 mm.*



*The **GK hektor** chair-lifter is the ideal extra for GK shower cubicles.*